
Bag of Tricks for Long-Tailed Multi-Label Classification on Chest X-Rays

Feng Hong*, Tianjie Dai*, Jiangchao Yao, Ya Zhang, Yanfeng Wang

Cooperative Medianet Innovation Center, Shanghai Jiao Tong University Shanghai AI Laboratory

Framework



- ► Image encoder, pretrained text encoder and query network.
 - The image encoder can be almost any well-known vision models. We preserve the visual features before pooling that can provide richer information in the subsequent query network.
 - During training, we freeze the text encoder to avoid overfitting.
 - After extracting the features of the two modalities, we use multiple transformer decoder layers for disease diagnosis, with label embedding as Query, and CXR image feature as Key and Value. The transformer decoder aims to combine textual information into corresponding image features, further improving predictive performance.

General Framework

 $S(x) = \Phi_{\text{query}}(\Phi_{\text{image}}(x), \Phi_{\text{text}}(\mathcal{Y}))$, where x denotes image, \mathcal{Y} denotes label set and $S(x) \in [0, 1]^{|\mathcal{Y}|}$ denotes the prediction of x.

LT-specific Designs

- ► Design 1: Textual Feature Extractor

 Two pretrained models: finetuned PubMedBERT or clinical-T5
- ► Design 2: Loss Function Reweighting

 Upweight loss of classes performed poorly on development set
- ► Design 3: Separate Classifier

For the last linear layer of the query network, we assign class-wise classifiers for each given class to preserve the independence of predictions between classes.

► Design 4: Data Augmentation

MixUp is a technique that randomly mixes up different images and labels.

► Design 5: Test-time Augmentation (TTA)

Multiple random transformations are applied on the test image and we test the model on all transformed images, and take the average as the output.

Design 6: Class-Wise Ensemble

For each class, we select best models based on development set and average their predictions.

► Design 7: Exogenous Data Replenishment

	Training Dataset	Amount	Label Space
	MIMIC-CXR	264849	26
	ChestXRay14	98637	14
	CheXpert	223414	14
•	Sum	586900	26
-			

Results

On development set

Table 1. Results of some different combinations of designs on the development set.									
Image Encoder	Text Encoder	Separate Classifier	Reweighting	MixUp	TTA	mAP			
ResNet-50	PubMedBERT					0.3187			
ResNet-50	PubMedBERT	\checkmark				0.3193			
ResNet-50	PubMedBERT	\checkmark			\checkmark	0.3220			
ResNet-50	PubMedBERT	\checkmark	\checkmark , 2		\checkmark	0.3273			
DenseNet-121	PubMedBERT	\checkmark	\checkmark , 2		\checkmark	0.3236			
ResNet-50	Clinical-T5	\checkmark	\checkmark , 3	\checkmark	\checkmark	0.3273			
ResNet-50	PubMedBERT	\checkmark	\checkmark , 3	\checkmark	\checkmark	0.3247			
ResNet-50	Clinical-T5	\checkmark	\checkmark , 2	\checkmark	\checkmark	0.3244			
ResNet-50	PubMedBERT	\checkmark	$\checkmark, 2$	\checkmark	\checkmark	0.3280			

- ► When separate classifier, reweighting(upweighting ratio=2), MixUp and TTA are simultaneously combined together with ResNet-50 and PubMedBERT, the best mAP score is achieved.
- ► The incorporation of a series of LT-specific designs help boost the chest X-rays disease diagnosis ability of the base general architecture from 0.3187 to 0.3280.

On testing set

Table 2. Results of model ensemble and data replenishment on the

test set.

test set.						
	Ensemble Method	Data Replenishment	mAP			
	Model-wise Ensemble		0.347			
	Class-wise Ensemble		0.347 0.348 0.349			
	Class-wise Ensemble	✓	0.349			

- ► Class-wise ensemble outperforms the model-wise ensemble, which directly averages the predictions of distinct models.
- ► External data replenishment further improves the performance, achieving the highest score among all attempts at a remarkable 0.349 mAP and ranking in the top five of ICCV CVAMD 2023 CXR-LT Competition.



